

HPP271 Hydrogen Peroxide Probe



				1	2	3	4	5
settings	Hydrogen Peroxide Probe		SPARE CODE	HPP271			O	O
	1 4..20 mA output signals for CH1	OFF			N			
		H2O2 0...2000 ppm			A			
		Dew Point Temperature (H2O and H2O2) +5...50 °C			J			
accessory	2 4..20 mA output signals for CH2	OFF			N			
		H2O2 0...2000 ppm			A			
		Dew Point Temperature (H2O and H2O2) +5...50 °C			J		O	
4 Cable	No Cable							N
	• <i>Probe Cables with Open Wires (1.5m)</i>		<i>254294SP</i>					<i>A</i>
	• <i>Probe Cables with Open Wires (3m)</i>		<i>254295SP</i>					<i>B</i>
	• <i>Probe Cables with Open Wires (5m)</i>		<i>254296SP</i>					<i>C</i>
	• <i>Probe Cables with Open Wires (10m)</i>		<i>254297SP</i>					<i>D</i>

NOTE: Selections in bold are included in the prices of the basic versions.
 • *Selections in italic are available at an extra price.*

LINKS: [User Guides](#)
[Datasheet](#)
[Vaisala.com product page](#)

Example of order code with typical settings:

		1	2	3	4	5
HPP271		A	N	O	B	O

Related Items

- Connection Cable M12-5F - Open Wires, Hygienic Design, 1.5m
- Connection Cable M12-5F - Open Wires, Hygienic Design, 3m
- Connection Cable M12-5F - Open Wires, Hygienic Design, 5m
- Connection Cable M12-5F - Open Wires, Hygienic Design, 10m
- Flat Cable, M12-5(M) - M12-5(F), 1m
- Filter Cap Porous PTFE 18,5mm HPP272 Spare Part
- Gland Set for Through-wall Installation, HPP271
- Flange for Through-wall Installation, HPP271
- Wall Mount for HPP271 and HPP272
- Probe Connection Cable, M12-5F - M12-5M, 10 m, Hygienic Design
- Probe Connection Cable, M12-5F - M12-5M, 1.5 m, Hygienic Design
- Probe Connection Cable, M12-5F - M12-5M, 3 m, Hygienic Design
- Probe Connection Cable, M12-5F - M12-5M, 5 m, Hygienic Design
- Indigo-compatible USB Adapter

Part number

- 254294SP
- 254295SP
- 254296SP
- 254297SP
- CBL210493SP
- DRW246363SP
- HPP271MOUNTINGSET1
- HPP271MOUNTINGSET2
- HPP272WALLMOUNT
- INDIGOCABLEHD10M
- INDIGOCABLEHD1M5
- INDIGOCABLEHD3M
- INDIGOCABLEHD5M
- USB2